## Generic Recon Chain and Toolboxes demo and inline debug

Hui Xue

NHLBI, NIH July 2020



#### In this session

- Motivation of generic chain
- MRD dimensions
- Triggering and on-the-fly recon
- Core data structures
- Walk through of key gadgets
- Interactive session set up Gadgetron in Ubuntu 20.04 from scratch and inline debug




#### Why the generic chain is needed

- To support common recon tasks
- To support deployment of Gadgetron
- Provide default implementation for reconstruction components
- Enable on-the-fly recon
- Demo and tutorial




RO, E1, CHA, PHS, SLC




RO, E1, CHA, SET, SLC, REP

T1



RO, E1, CHA, SET, SLC



RO, E1, E2, CHA, SLC





#### MRD data dimensions

We don't want to write a separate chain for every imaging app

MR kspace data acquisition is sequential in general (e.g. compared to 3D ultrasound)

MRD data format encodes this data acquisition scheme

```
uint16 t version;
                                                                                                                                                               /**< First unsigned int indicates the version */
                                                                                                             uint16_t data_type;
                                                                                                                                                               /**< e.g. unsigned short, float, complex float, etc. */
                                                                                                             uint64 t flags;
                                                                                                                                                               /**< bit field with flags */
                                                                                                             uint32_t measurement_uid;
                                                                                                                                                               /**< Unique ID for the measurement */
                                                                                                             uint16_t matrix_size[3];
                                                                                                                                                               /**< Pixels in the 3 spatial dimensions */
                                                                                                             float field_of_view[3];
                                                                                                                                                               /**< Size (in mm) of the 3 spatial dimensions */
                                                                                                                                                               /**< Number of receive channels */
typedef struct ISMRMRD EncodingCounters {
                                                                                                             uint16_t channels;
                                                                                                             float position[3]:
                                                                                                                                                               /**< Three-dimensional spatial offsets from isocenter */
    uint16 t kspace encode step 1;
                                              /**< e.g. phase encoding line number */</pre>
                                                                                                             float read dir[3];
                                                                                                                                                               /**< Directional cosines of the readout/frequency encoding */
    uint16 t kspace encode step 2;
                                             /**< e.g. partition encoding number */
                                                                                                             float phase_dir[3];
                                                                                                                                                               /**< Directional cosines of the phase */
                                              /**< e.g. signal average number */
    uint16_t average;
                                                                                                                                                               /**< Directional cosines of the slice direction */
                                                                                                             float slice_dir[3];
                                                                                                             float patient_table_position[3];
                                                                                                                                                               /**< Patient table off-center */
                                              /**< e.g. imaging slice number */</pre>
    uint16_t slice;
                                                                                                             uint16 t average;
                                                                                                                                                               /**< e.g. signal average number */
    uint16 t contrast;
                                              /**< e.g. echo number in multi-echo */
                                                                                                             uint16_t slice;
                                                                                                                                                               /**< e.g. imaging slice number */
    uint16_t phase;
                                              /**< e.g. cardiac phase number */</pre>
                                                                                                             uint16 t contrast:
                                                                                                                                                               /**< e.g. echo number in multi-echo */
                                                                                                             uint16 t phase;
                                                                                                                                                               /**< e.g. cardiac phase number */
    uint16 t repetition;
                                              /**< e.g. dynamic number for dynamic scanning */
                                                                                                             uint16_t repetition;
                                                                                                                                                               /**< e.g. dynamic number for dynamic scanning */
    uint16 t set;
                                              /**< e.g. flow encoding set */
                                                                                                             uint16_t set;
                                                                                                                                                               /**< e.g. flow encodning set */
                                              /**< e.g. segment number for segmented acquisitio
    uint16 t segment;
                                                                                                             uint32_t acquisition_time_stamp;
                                                                                                                                                               /**< Acquisition clock */
    uint16 t user[ISMRMRD USER INTS]; /**< Free user parameters */</pre>
                                                                                                             uint32_t physiology_time_stamp[ISMRMRD_PHYS_STAMPS]; /**< Physiology time stamps, e.g. ecg, breathing, etc. */
                                                                                                             uint16_t image_type;
                                                                                                                                                               /**< e.g. magnitude, phase, complex, real, imag, etc. */
} ISMRMRD EncodingCounters;
                                                                                                             uint16_t image_index;
                                                                                                                                                               /**< e.g. image number in series of images */
                                                                                                             uint16_t image_series_index;
                                                                                                                                                               /**< e.g. series number */
                                                                                                             int32_t user_int[ISMRMRD_USER_INTS];
                                                                                                                                                               /**< Free user parameters */
                                                                                                             float user_float[ISMRMRD_USER_FLOATS];
                                                                                                                                                               /**< Free user parameters */
                                                                                                             uint32_t attribute_string_len;
                                                                                                                                                               /**< Length of attributes string */
                                                                                                          1 TOURING Tournelles dans
```

typedef struct ISMRMRD\_ImageHeader {



#### Order of acquisition to fill multidimensional array does matter

E.g. Cine imaging – RO, E1, CHA, PHS, SLC


Let's look at an example ...


meas\_MID838\_PK\_rt\_test\_2slice\_FID22519



#### On-the-fly recon

- No need to wait, but start when data was available
- Should be flexible








#### Consider the inline behavior


On-the-fly recon is popular among MR technicians. Waiting time: 2s with on-the-fly, 20s without on-the-fly



National Heart, Lung, and Blood Institute

#### Generic recon chain

#### Key idea/features:





#### Generic recon chain: core data structure

#### IsmrmrdAcquisitionBucket

gadgetron\toolboxes\mri\_core\mri\_core\_acquisition\_bucket.h

 Un-ordered buffer of incoming kspace lines, subject to sorting\_dimension

```
·class·IsmrmrdAcquisitionData
· · · · GadgetContainerMessage < ISMRMRD:: AcquisitionHeader > * · head ;

· h · GadgetContainerMessage( · hoNDArray( · std::complex(float) · > · > * · data;

· GadgetContainerMessage(·hoNDArray(·float·>·>·*·traj;
· };
 ·class · IsmrmrdAcquisitionBucketStats
 ٠ {
 · · · public:
 · · · · · // · Set · of · labels · found · in · the · data · or · ref · part · of · a · bucket
 · · · · · //11D, · fixed · order · [RO, · E1, · E2, · CHA, · SLC, · PHS, · CON, · REP, · SET, · SEG, · AVE]
 ....std::set<uint16 t>.kspace encode step 1;
 ....std::set<uint16_t>.kspace_encode_step_2;
 ....std::set<uint16 t>.slice;
 ....std::set<uint16 t>.phase;
 ....std::set<uint16 t>.contrast;
 ....std::set<uint16_t>.repetition;
 ....std::set<uint16 t>.set;
 ....std::set<uint16 t>.segment;
 ....std::set<uint16 t>.average;
 · };
```



#### Generic recon chain: core data structure

#### IsmrmrdReconData

An IsmrmrdAcquisitionBucket is converted to one or more IsmrmrdReconData seven dimensional array: [RO, E1, E2, CHA, N, S, SLC]

·struct · IsmrmrdReconBit

```
·struct·IsmrmrdReconData

·{

·public:

·rstd::vector<IsmrmrdReconBit>·rbit_;

·};
```

- Every encoding space corresponds to a ReconBit
- A ReconBit has data and ref
- Data and ref has 7D arrays for kspace/traj/density/headers...

```
public:
  ·IsmrmrdDataBuffered·data ;
  ·boost::optional<IsmrmrdDataBuffered>·ref ;
·struct · IsmrmrdDataBuffered
·public:
·/·//TD, ·fixed·order·[E0, ·E1, ·E2, ·CHA, ·N, ·S, ·LOC]
...hoNDArray<.std::complex<float>.>.data ;
···//7D, ·fixed·order·[TRAJ, ·E0, ·E1, ·E2, ·N, ·S, ·LOC]
boost::optional<hoNDArray<float>> trajectory_;
·/·/·6D, ·density·weights·[E0, ·E1, ·E2, ·N, ·S, ·LOC]
. boost::optional<hoNDArray<float>.>.density_;
...//5D, fixed order [E1, E2, N, S, LOC]
...hoNDArray<.ISMRMRD::AcquisitionHeader.>.headers ;
SamplingDescription sampling;
```

#### Generic recon chain: core data structure

IsmrmrdImageArray data, header, meta fields Each [RO, E1, E2, CHA] shares one header



#### Key gadgets: Triggering

```
·class·EXPORTGADGETSMRICORE·AcquisitionAccumulateTriggerGadget::
·public·Gadgetron::Gadget10f2<ISMRMRD::AcquisitionHeader, ·ISMRMRD::ISMRMRD_WaveformHeader.>
. . . {
· · · public:
· · · · · · GADGET_DECLARE(AcquisitionAccumulateTriggerGadget);
····typedef·std::map<·unsigned·short·int, GadgetContainerMessage<IsmrmrdAcquisitionBucket>*·>·map_type_;
····virtual·~AcquisitionAccumulateTriggerGadget();
...int·close(unsigned·long·flags);
···protected:
·····GADGET_PROPERTY_LIMITS(trigger_dimension, ·std::string, ·"Dimension · to ·trigger · on", ·"",

→ → ·····GadgetPropertyLimitsEnumeration,

→ 
→ 
·····"kspace_encode_step_1",
      "kspace_encode_step_2",
      → ·····"average",
       → "contrast",
       → ·····"phase",
       "repetition",
       → ·····"set",
       → ·····"segment",
       → "user_0",
       ....."user_1",
.....n_acquisitions",
   → → ·····");
····GADGET_PROPERTY_LIMITS(sorting_dimension, std::string, "Dimension to sort by", ."", .
   → GadgetPropertyLimitsEnumeration,
     ....."kspace_encode_step_1",

> > ..... "kspace_encode_step_2",
→ → ·····"average".
```

```
<!---Data accumulation and trigger gadget -->

<gadget>

....<name>AccTrig</name>

....<dll>gadgetron_mricore</dll>

....<classname>AcquisitionAccumulateTriggerGadget</classname>

....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
....
```

#### gadgets\mri\_core\ AcquisitionAccumulateTriggerGadget

Input	Kspace line as ISMRMRD::AcquisitionHeader and hoNDArray
Output	IsmrmrdAcquisitionBucket
Key paras	trigger_dimension : once this dimension changes, send the buffered data to next gadget
Notes	<ul> <li>If trigger_dimension is not set, trigger only once in close() to send all data to next</li> <li>sorting_dimension, if set, sort the buffered kspace line before sending</li> </ul>



#### Key gadgets: BucketToBuffer

```
·class·EXPORTGADGETSMRICORE·BucketToBufferGadget·:·public·Gadget1<IsmrmrdAcquisitionBucket>
···public:
·····GADGET DECLARE(BucketToBufferGadget);
·····BucketToBufferGadget();
····virtual ~BucketToBufferGadget();
....int.close(unsigned.long.flags);
GADGET PROPERTY LIMITS(N dimension, std::string, "N-Dimensions", "",
······GadgetPropertyLimitsEnumeration,
...."average",
.....contrast",
    ·····phase",
repetition",
.....segment",
.....slice",
.,....);
·····GADGET_PROPERTY_LIMITS(S_dimension, ·std::string, ·"S-Dimensions", ·"", ·
GadgetPropertyLimitsEnumeration,
....."average",
     ·····"contrast",
· repetition",
.....segment"
....."slice",
.....;
·····GADGET PROPERTY(split slices, bool, "Split slices", false);
....GADGET PROPERTY(ignore segment, bool, "Ignore segment", false);
·····GADGET_PROPERTY(verbose, bool, "Whether to print more information", false);
```

Bucket has AVE=3, N is Phase and S is SET


Then three ReconData will be sent out for ave=0,1,2

· <gadget></gadget>	
···· <name< th=""><th>e&gt;BucketToBuffer</th></name<>	e>BucketToBuffer
···· <dl1></dl1>	gadgetron_mricore
····· <clas< td=""><td>ssname&gt;BucketToBufferGadget</td></clas<>	ssname>BucketToBufferGadget
···· <prop< td=""><td>perty&gt;<name>N_dimension</name><value>contrast</value></td></prop<>	perty> <name>N_dimension</name> <value>contrast</value>
···· <prop< td=""><td>perty&gt;<name>S_dimension</name><value>average</value></td></prop<>	perty> <name>S_dimension</name> <value>average</value>
····· <prop< td=""><td>perty&gt;<name>split_slices</name><value>false</value></td></prop<>	perty> <name>split_slices</name> <value>false</value>
····· <prop< td=""><td>perty&gt;<name>ignore_segment</name><value>true</value></td></prop<>	perty> <name>ignore_segment</name> <value>true</value>
···· <prop< td=""><td>perty&gt;<name>verbose</name><value>true</value></td></prop<>	perty> <name>verbose</name> <value>true</value>
·	•

gadgets\mri_core\
BucketToBufferGadget

Input	IsmrmrdAcquisitionBucket
Output	One or more IsmrmrdReconData
Key paras	N_dimension, S_dimension: [RO, E1, E2, CHA, N, S, SLC]
Notes	<ul> <li>Every incoming bucket is split into one or more ReconData</li> <li>Separated into imaging and ref arrays</li> <li>Optionally, include trajectory, density func</li> <li>Include acquisition headers</li> </ul>

#### Key gadgets: Prepare reference



#### Key gadgets: Coil compression

<pre>class-EXPORTGADGETSMRICORE-GenericReconEigenChannelGadget::- {     public:</pre>	oublic-GenericReconDataBase	<pre><!--Coil.compression--></pre>	
//////-parameters-to-control-the-reconstruction//////-compute-KLT-coefficients			
///whether.to:average.all.N.for.coefficient.computatio/// for.the.interleaved.mode.the.sampling.times.will.bGADGET_PROPERTY(average_all.pef_N,.bool,."Whether.to:av/// whether.to:average.all.S.for.coefficient.computatioGADGET_PROPERTY(average_all_ref_S,.bool,."Whether.to:av///if.update_eigen_channel_coefficients==true,.every.i/// and.the.older.one.will.be.replaced	Input	IsmrmrdReconData	
///if·update_eigen_channel_coefficients==false, the KLGADGET_PROPERTY(update_eigen_channel_coefficients, bool///optionally, upstream-coil-compression-can-be-applie//if·upstream_coil_compression==true, only kept-chann///no.matter.whether.upstream_coil_compression-is-trueGADGET_PROPERTY(upstream_coil_compression, bool, "Wheth	Output	IsmrmrdReconData	
//·the·logic·here·is·that·if·upstream_coil_compression///if·upstream_coil_compression_num_modesKept<=0·and·ukeep//·the·first·N·and·first·S·will·be·used·to·compute·numGADGET_PROPERTY(upstream_coil_compression_thres,.doubleGADGET_PROPERTY(upstream_coil_compression_num_modesKept	Key paras	Control how to compute coil compression coefficients	
	Notes	<ul> <li>Compute coefficients on ref data</li> <li>Control how many channels to keep with thres and num_modesKept</li> <li>Different coefficients for different SLC</li> <li>Output kspace are always in eigen channel (high energy first, along CHA)</li> </ul>	

and Blood Institute

#### Key gadgets: Recon

class·EXPORTGADGETSMRICORE·GenericRecond { public: GADGET_DECLARE(GenericReconCartesian typedef·GenericReconGadget·BaseClass typedef·Gadgetron::GenericReconCarte	nGrappaGadget);	<pre><dll>gadgetron_mricore</dll><classname>GenericReconCartesianGrappaGadget</classname><!--image·series--><pre>property&gt;<name>image_series</name><value>@</value><!--Coil·map·estimation, Inati-or-Inati_Iter--><pre>property&gt;<name>coil_map_algorithm</name><value>Inati</value></pre>/property&gt;</pre></pre>
GenericReconCartesianGrappaGadget(  GenericReconCartesianGrappaGadget  ///  ///  parameters to control the reco  ///  ///  GADGET_PROPERTY(send_out_sfactor,  GADGET_PROPERTY(send_out_snr_map,	Gen	gadgets\mri_core\ GenericReconGadget, enericReconCartesianGrappaGadget, ericReconCartesianGrappaAlGadget, enericReconCartesianGrappaAlGadget. enericReconCartesianSpiritGadget
///·Grappa·parametersGADGET_PROPERTY(grappa_kSize_RO, i	Input	IsmrmrdReconData
GADGET_PROPERTY(grappa_kSize_E2, i GADGET_PROPERTY(grappa_reg_lamda, GADGET_PROPERTY(grappa_calib_over_	Output	IsmrmrdImageArray
	Key paras	Parameters to control recon, e.g. grappa, kernel size, regularization strength
16	Notes	<ul> <li>Key task is to convert ReconData to array of images</li> <li>Follow 7D convention</li> <li>User should implement their recon methods here</li> </ul>

## Key gadgets: Partial Fourier handling, Kspace filter

**Notes** 

· Kspace·fi</th <th>Itering&gt;</th>	Itering>
- <gadget></gadget>	
···· <name>KSpac</name>	eFilter
····· <dll>gadget</dll>	tron_mricore
····· <classname< td=""><td>GenericReconKSpaceFilteringGadget</td></classname<>	GenericReconKSpaceFilteringGadget
···· parame</td <td>eters for debug and timing&gt;</td>	eters for debug and timing>
···· <property></property>	<pre>(name&gt;debug_folder&lt;<value></value></pre>
····· <property></property>	<pre>(name&gt;perform_timing<value>false</value></pre>
···· <property></property>	<pre>cname&gt;verbose<value>false</value></pre>
···· if-in</td <td>coming images have this meta field, it will not be processed&gt;</td>	coming images have this meta field, it will not be processed>
····· <property></property>	<pre>cname&gt;skip_processing_meta_field<value>Skip_processing_after_recon</value></pre>
···· parame</td <td>eters for kspace filtering&gt;</td>	eters for kspace filtering>
····· <property></property>	<pre>(name&gt;filterRO<value>Gaussian</value></pre>
···· <property></property>	<pre>(name&gt;filterRO_sigma<value>1.0</value></pre>
···· <property></property>	<pre>cname&gt;filterRO_width<value>0.15</value></pre>
···· <pre>property&gt;</pre>	<pre>cname&gt;filterE1<value>Gaussian</value></pre>
···· <pre>property&gt;</pre>	<pre>(name&gt;filterE1_sigma<value>1.0</value></pre>
···· <property></property>	<pre>cname&gt;filterE1_width<value>0.15</value></pre>
···· <pre>property&gt;</pre>	<pre>(name&gt;filterE2<value>Gaussian</value></pre>
····· <property></property>	<pre>(name&gt;filterE2_sigma<value>1.0</value></pre>
···· <property></property>	<pre>(name&gt;filterE2_width<value>0.15</value></pre>
·	


## gadgets\mri\_core\ GenericReconPartialFourierHandlingGadget GenericReconPartialFourierHandlingFilterGadget GenericReconPartialFourierHandlingPOCSGadget Input IsmrmrdImageArray Output IsmrmrdImageArray Key paras Decide which gadget to use, filter strength, POCS

Support 2D and 3D

Detect Partial Fourier from encoding space parameters

gadgets\mri_core\ GenericReconKSpaceFilteringGadget		
Input	IsmrmrdImageArray	
Output	IsmrmrdImageArray	
Key paras	filterRO, filterE2, filterE3	
Notes	<ul> <li>Separable filter for RO/E1/E2</li> <li>Default Gaussian filter</li> <li>Good trade-off of main lobe width and remote lobe suppression</li> </ul>	
	and Blood Institute	


#### Kspace filter



Slides from Dr. Peter Kellman, NHLBI, NIH



#### Kspace filter: Comparison



Gaussian, sigma=1.5pixel

Turkey, weak on the scanner

Turkey, medium on the scanner

Slides from Dr. Peter Kellman, NHLBI, NIH



#### Key gadgets: FOV adjust and scaling

# GenericReconFieldOfViewAdjustmentGadget Input IsmrmrdImageArray Output IsmrmrdImageArray Key paras Ismrmrd protocol Notes Support 2D and 3D Make sure output images are in correct size/FOV

<pre>class-EXPORTGADGETSMRICORE-GenericReconImageArrayScalingGadget-:-public-GenericReconImageBase</pre>
<pre>typedef·float·real_value_type;typedef·std::complex<real_value_type>-ValueType;typedef·ValueType·T;</real_value_type></pre>
<pre>typedef-GenericReconImageBase-BaseClass;GenericReconImageArrayScalingGadget();~GenericReconImageArrayScalingGadget();</pre>
<pre>'''' '''' parameters to control the reconstruction '''' '''' image scaling 'GADGET PROPERTY(use_constant_scalingFactor, bool, "Whether to use constrant scaling; if not, the auto-scaling computed only ONCE", true); 'GADGET PROPERTY(scalingFactor, float, "Default scaling ratio", 4.0); 'GADGET PROPERTY(min_intensity_value, int, "Minimal intensity value for auto image scaling", 64); 'GADGET PROPERTY(max_intensity_value, int, "Maximal intensity value for auto image scaling", 4095); 'GADGET PROPERTY(auto_scaling_only_once, bool, "Whether to compute auto-scaling factor only once; if false factor is computed for every incoming image array", true);</pre>

	gadgets\mri_core\ GenericReconlmageArrayScalingGadget		
	Input	IsmrmrdImageArray	
	Output	IsmrmrdImageArray	
-5(	Key paras	use_constant_scalingFactor, auto_scaling_only_once	
al:	Notes	Recommendation  • Use constant scaling if SNR unit recon  • Only scale once for multiple image array in one scan	

and Blood Institute

#### **Example: RealTime Cine**


#### gadgets/mri\_core/config/Generic\_Cartesian\_Grappa\_RealTimeCine.xml

```
.<!--Data-accumulation-and-trigger-gadget--->
<gadget>
····<name>AccTrig</name>
·····<dll>gadgetron mricore</dll>
·····<classname>AcquisitionAccumulateTriggerGadget</classname>
.....property><name>trigger dimension</name><value>slice</value></property>
.....property><name>sorting dimension</name></value>
.</gadget>
<gadget>
····<name>BucketToBuffer</name>
·····<dll>gadgetron mricore</dll>
·····<classname>BucketToBufferGadget</classname>
.....property><name>N_dimension</name><value>phase</value>

....property><name>S dimension</name><value>set</value>
/property>
.....property><name>split_slices</name><value>true</value>
·····property><name>ignore_segment</name><value>true</value></property>
</gadget>
·<!-- Prep ref ·-->
····<name>PrepRef</name>
·····<dll>gadgetron mricore</dll>

....<classname>GenericReconCartesianReferencePrepGadget</classname>
····<!---parameters.for.debug.and.timing.-->
.....property><name>debug_folder</name><value></value></property>
....property><name>perform timing</name><value>true</value>
.....continue of the continue of the
····<!---averaging across repetition --->
·····property><name>average_all_ref_N</name><value>true</value></property>
····<!--every.set.has.its.own.kernels.-->
.....property><name>average_all_ref_S</name><value>false</value>
····<!--·whether always to prepare ref if no acceleration is used -->
·····property><name>prepare_ref_always</name><value>true</value></property>
</gadget>
```

```
·<!--·Coil·compression·-->
<gadget>
····<name>CoilCompression</name>
·····<dll>gadgetron_mricore</dll>
·····<classname>GenericReconEigenChannelGadget</classname>
····<!---parameters.for.debug.and.timing.-->
·····<property><name>debug_folder</name><value></value></property>
·····property><name>perform timing</name><value>true</value>
·····property><name>verbose</name><value>true</value></property>
·····property><name>average_all_ref_N</name><value>true</value>
·····property><name>average_all_ref_S</name><value>true</value>
····<!---Up·stream·coil·compression·-->
·····property><name>upstream_coil_compression</name><value>true</value>
·····roperty><name>upstream_coil_compression_thres</name><value>-1</value>
 ·····continued in the continue of the con
·<!---Recon--->
· · · · <name>Recon</name>
·····(dll)gadgetron mricore(/dll)
····<classname>GenericReconCartesianGrappaGadget</classname>
····<!---image-series--->
·····property><name>image_series</name><value>0</value>
····<!---Coil·map·estimation, ·Inati·or·Inati_Iter--->
·····<property><name>coil_map_algorithm</name><value>Inati</value></property>
····<!--Down-stream-coil-compression--->
·····<property><name>downstream_coil_compression</name><value>true</value></property>
·····<property><name>downstream_coil_compression_thres</name><value>0.002</value></property>
·····<property><name>downstream_coil_compression_num_modesKept</name><value>0</value></property>
····<!---parameters.for.debug.and.timing.-->
·····property><name>debug_folder</name><value></value></property>
·····perty><name>perform_timing</name><value>true</value>
·····property><name>verbose</name><value>true</value></property>
 ····<!---whether-to-send-out-gfactor--->
·····property><name>send_out_gfactor</name><value>true</value></property>
 </gadget>
·<!-- Partial fourier handling ·-->
····<name>PartialFourierHandling</name>
····<dll>gadgetron_mricore</dll>
```



### Demo session Generic chain in ubuntu

- Let's start from scratch
  - Install dependencies
  - Clone Gadgetron, ismrmrd etc.
  - cmake/ccmake
  - Compile debug version and install
  - Inline debug/code walk through with Eclipse
  - Walk through key gadgets in generic chain

For step-by-step instructions:

https://github.com/gadgetron/gadgetron/wiki/Visual-debug-Gadgetron-in-Ubuntu-using-Eclipse





